


BD Venflon™ Pro Safety IV cannula

Safety is in sight.

Ported safety peripheral IV catheter.



1 Designed for smooth insertion.
A back-cut needle, smooth catheter transition tapered tip and proprietary lube system designed to help minimise insertion discomfort.*¹

2 Safeguard veins with longer dwell times.
Proprietary catheter material, BD Vialon™, supported by multiple clinical studies.^{2,3} Its unique properties allow to soften up to 70% in the vein¹, enabling longer dwell times and reducing the chance of mechanical phlebitis by up to 50%.²

3 Comfortable securement.
Flexible wings stabilise the cannula within the vein and provide comfortable securement for the patient.

4 Familiar design.
Familiar ported design with ‘snap’ cap aids non-touch technique and helps prevent the port from being accidentally opened.

5 Care for the caregiver.
The BD Venflon™ Pro Safety safety IV catheter utilises a passive safety system to enclose the needle, providing protection from needlestick injury and resulting exposure.⁴

Captures residual blood safely within the device, minimising the risk of unanticipated blood spatter during and after activation.⁵

Power injectable catheter.

The BD Venflon™ Pro Safety catheters are power injectable to a maximum of 325 PSI.⁶

1. C-1974 C-0030 Infusion Claim 4 - BD Vialon biomaterial (v1.0). 05/2020.

2. Maki DG, Ringer M. Risk factors for infusion-related phlebitis with small peripheral venous catheters. *Annals of Internal Medicine*. 1991;114:845-854.

3. Stanley MD, Meister E, Fuschuber K. Infiltration during intravenous therapy in neonates: comparison of Teflon® and Vialon® catheters. *Southern Medical Journal*. 1992;85:883-886.

4. TECHNICAL DATA SHEET BD Venflon™ Pro Safety IV catheter. Revision May 2018.

5. Rosenthal VD and Hughes G. Fluid dispersal from safety cannulas: An in vitro comparative test. *Am J Infect Control* 2015.BD-8674

6. BD Venflon & Venflon Pro Safety with BD Instaflash Power Injection Flow Rates White Paper

* Compared to a non-tapered tip

BD Switzerland Särl, Terre Bonne Park – A4 Route de Crassier 17, 1262 Eysins, Switzerland

bd.com